A560/560PC A580/580PC A590/590PC SPECTROPHOTOMETER

Features:

- Large LCD screen(480*800) with 65 thousand true color and support U disk data saving
- Imported & Environmental Deuterium lamp, preventing from Ozone inhalation
- Data output: Data can be exported to U disk (USB Output)
- Print port: Instrument can connect PCL printer for printing A4 paper (USB Port)
- O Adopted high-class grating with wholly hermetic light path design, to ensure the instrument has the super low stray light
- Upgrade package can be downloaded from AOE web with S/N, support U disk one-key upgrade
- Real-time monitoring the lifetime of Deuterium lamp and Tungsten lamp with advanced system
- Pre-aligned design ensures the user can change lamps conveniently
- With GLP self-check function, check the wavelength accuracy and Photometry accuracy, can provide test report and power-off protection
- Wavelength calibration, wavelength setting, change lamp source and dark current calibration automatically
- With SiO₂ coating optical mirror, reducing the pollution from outside fully

Basic Functions:

1. Photometry (Basic)

Test Abs., Transmittance and Energy by fixed wavelength

2. Quantitative (Qty)

Linear fit and Linear fit through zero two modes

- a. Coefficient, Standard Sample input and Standard Sample read three modes to establish standard curve
- b. Establish A=K1*C+K0, can search original data, graph curve, parameters settings
 c. Can save 240 group curves, can test 240 data in each curve

3. Kinetics (Kins)

Used for time course scanning or reaction rate calculations $\triangle A/t$, can search all data

4. Multi-Wavelength (MultiWL)

Can test Transmittance and Abs. with 8 different wavelengths at most

5. Spectrum (Spec)

User can set the scan range from 190nm to 1100nm to test the max. Abs. peak value, can do derivation, arithmetical operations to the graph.



FAITHFUL

d. Double wavelength, Triple wavelength test functions.

Specifications

6. BIO (DNA/Protein)

ONA/Protein. (Optional:UV, Lowry, BCA, CBB and Biuret)

Model	A560 A560PC	A580 A580PC	A590 A590PC
Optical System	Double Beam, Grating 1200 lines/mm		
Wavelength Range	190–1100nm		
Wavelength Range	1.8nm	1nm	0.5. 1. 2. 4, 5nm
Wavelength Accuracy	±0.3nm		
Wavelength Repeatability	≤ 0.2nm		
Photometric Accuracy	$\pm 0.002A (Q-0.5Abs) \pm 0.004A (0.5=1.0Abs) \pm 0.3\% T (0-100\% T)$		
Photometric Repeatability	$0.001Abs (0-0.5Abs), 0.002Abs (0.5-1.0Abs), \le 0.1\% T (0-100\% T)$		
Stray Light	≤0.03% T @ 360nm;220nm		
Stability	± 0.0003A / h @ 500nm		
Baseline Flatness	± 0.0005A		
Noise	0.0002Abs		
Display	65 thousand ture color 7 inch TFT LCD(480*800)		
Photometric Mode	T, A, C, E		
Photometric Range	0-200% T,-0.301-3.0A		
Dete tor	Silicon Photodiode		
Light Source	Deuterium Lamp, Tungsten Lamp		
Input	Membrare Keypad		
Output	USB-A* 2 Print and data output: USB-B Connect PC		
Compartment	Optional 8 Auto Cell holder, Solid Sample Holder, Micro Cell Holder, 10-100mm Cell Holder		