

# FAITHFUL® FAITHFUL INSTRUMENT (Hebei) CO.,Ltd

#### CLARIFY TESTER CM SERIES

They test clarify of injection and bottled medicine liquid.

### Applicable standards

- National standard(Chinese codex )/2015
- Health department standard

Corporate standard(Clarify tester)

## FCM-1 CLARIFY TESTER

#### Features

- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- Flash e liminated. Illumination is adjustable.
- Time and alarm automatically. Test continuously.
- Used by one person.

## Specifications

- Lamp three-primary-color fluorescent lamp
- Illumination range (1000~4000) LUX
- Power 220V / 50Hz/ 30W
- $(1\sim 99) s$ ●Time range



Dimension (66\*13\*50) cm

## FCM-2 CLARIFY TESTER



#### Features

- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.
- Flash eliminated. Illumination is adjustable.
- Time and alarm automatically. Test continuously.
- Used by two person at the same time.

## Specifications

- Lamp three-primary-color fluorescent lamp
- Illumination range (1000~4000) LUX
- Time range  $(1\sim 99) s$
- Power 220V / 50Hz/ 50W
- Dimension (66\*38\*50) cm